

NOVO Respiratory Assist Device (RAD)



Positive Pressure and Oxygen Therapy System with
Protection for Caregivers

WIP Snapshot 04/19/2020

NOVO ENGINEERING®
Better products to market. *Faster.*

NOTICE: This device is under development. It has not been tested on humans or submitted for regulatory clearance/approval

PROJECT OBJECTIVES

- Provide practical **oxygen-enriched, positive airway pressure** therapy to patients with respiratory insufficiency at home or in a step-down facility.
- Incorporate safety features that **protect health care providers from aerosolized virus**



NOTE: We blocked the vents on this vented mask for testing purposes

PRESUMED BENEFITS TO PATIENTS AND CAREGIVERS



○ PATIENTS:

- Non-invasive single or bi-level positive pressure therapy
- Oxygen enrichment
- Delay or avoid progressive hypoxemia and invasive ventilation
- Accelerate weaning from ventilator
- Minimizes backpressure during cough or sneeze
- Constant flow of fresh air and oxygen purges CO₂
- **Designed for home use or step-down facilities**



○ CAREGIVERS:

- **Prevents cross infection by containing and filtering all exhaled breath through viral filters**
- Works with standard CPAP, BPAP, or VPAP equipment
- Adaptable to any full-face, **non-vented mask**, with nonelastic headgear
- Mask pressure limiting valve prevents mask from leaking during cough or sneeze
- Can be combined with nebulizers, suction catheters, and other respiratory care devices
- Simple setup and operation with independent control over PEEP pressure and flow rate
- Safe, convenient method to drain condensate from the system

ADDITIONAL FEATURES IN PROCESS




- ◉ Improved CO2 purging
- ◉ Integration of a suction catheter
- ◉ Anti-asphyxiation valve
- ◉ Optimization of the over-pressure valve
- ◉ Adjustable PEEP and pressure relief settings
- ◉ Injection moldable version of the design



INITIAL SYSTEM TESTING




NOVO Respiratory Assist Device (cm of H₂O)
2020 04 02 Pressure Data










RESPIRATORY ASSIST WITH OXYGEN, EXHALE FILTER, COUGH VALVE



BREATHING CYCLE



FLOW SIMULATIONS



- Use 3D-printed, screw machined, and die-cut parts for early versions
- Injection molding and die cut processes for higher-volume production
- For purchased items, specify standardized components available from multiple suppliers
- Recruit (medical) manufacturing operations to build with quality control
- Pair manufacturers with hospital systems close by for prescription & patient monitoring



WHAT WE NEED TO CONTINUE THIS WORK

- Clinicians to use the system under an informed consent agreement and provide feedback/support
- **Funding** to continue the development of the second-generation design.



THANK YOU

For more information contact:

CovidResponse@novoengineering.com

Dan Kline: dank@novoengineering.com

NOVO ENGINEERING®